



MAGEA5 Antibody, HRP conjugated

Product Code	CSB-PA013331LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43359
Immunogen	Recombinant Human Melanoma-associated antigen 5 protein (1-124AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Melanoma-associated antigen 5 (Cancer/testis antigen 1.5) (CT1.5) (MAGE-5 antigen), MAGEA5, MAGE5
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	MAGEA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.