



MAP3K12 Antibody, FITC conjugated

Product Code	CSB-PA013420LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12852
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 12 protein (101-400AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 12 (EC 2.7.11.25) (Dual leucine zipper bearing kinase) (DLK) (Leucine-zipper protein kinase) (ZPK) (MAPK-upstream kinase) (MUK) (Mixed lineage kinase), MAP3K12, ZPK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP3K12