



MAP3K12 Antibody, Biotin conjugated

Product Code	CSB-PA013420LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12852
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 12 protein (101-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 12 (EC 2.7.11.25) (Dual leucine zipper bearing kinase) (DLK) (Leucine-zipper protein kinase) (ZPK) (MAPK-upstream kinase) (MUK) (Mixed lineage kinase), MAP3K12, ZPK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP3K12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.