




MAP3K7 Antibody

Product Code	CSB-PA013429LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43318
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 7 protein (101-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 7 (EC 2.7.11.25) (Transforming growth factor-beta-activated kinase 1) (TGF-beta-activated kinase 1), MAP3K7, TAK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP3K7

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA013429LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.