



MPDU1 Antibody, Biotin conjugated

Product Code	CSB-PA014744LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75352
Immunogen	Recombinant Human Mannose-P-dolichol utilization defect 1 protein (1-36AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mannose-P-dolichol utilization defect 1 protein (Suppressor of Lec15 and Lec35 glycosylation mutation homolog) (SL15), MPDU1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MPDU1