



MPDU1 Antibody, Biotin conjugated

Product Code	CSB-PA014744LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75352
Immunogen	Recombinant Human Mannose-P-dolichol utilization defect 1 protein (1-36AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mannose-P-dolichol utilization defect 1 protein (Suppressor of Lec15 and Lec35 glycosylation mutation homolog) (SL15), MPDU1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MPDU1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.