



MPZL2 Antibody, FITC conjugated

Product Code	CSB-PA014776LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60487
Immunogen	Recombinant Human Myelin protein zero-like protein 2 protein (30-140AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myelin protein zero-like protein 2 (Epithelial V-like antigen 1), MPZL2, EVA EVA1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MPZL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.