



MRPL14 Antibody, Biotin conjugated

Product Code	CSB-PA014819LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P1L8
Immunogen	Recombinant Human 39S ribosomal protein L14, mitochondrial protein (1-145AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	39S ribosomal protein L14, mitochondrial (L14mt) (MRP-L14) (39S ribosomal protein L32, mitochondrial) (L32mt) (MRP-L32) (Mitochondrial large ribosomal subunit protein uL14m), MRPL14, MRPL32 RPML32
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MRPL14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.