



MTMR11 Antibody, HRP conjugated

Product Code	CSB-PA015176LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A4FU01
Immunogen	Recombinant Human Myotubularin-related protein 11 protein (1-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myotubularin-related protein 11 (Cisplatin resistance-associated protein) (hCRA), MTMR11
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MTMR11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.