



# MUC2 Antibody, Biotin conjugated

Product Code	CSB-PA015222LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02817
Immunogen	Recombinant Human Mucin-2 protein (36-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mucin-2 (MUC-2) (Intestinal mucin-2), MUC2, SMUC
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MUC2