




NADK Antibody

Product Code	CSB-PA015409LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95544
Immunogen	Recombinant Human NAD kinase protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NAD kinase (EC 2.7.1.23) (Poly(P)/ATP NAD kinase), NADK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NADK

Image



Western Blot

Positive WB detected in: SH-SY5Y whole cell lysate

All lanes: NADK antibody at 1.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 50, 64, 47 kDa

Observed band size: 47 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.