



NARS Antibody, Biotin conjugated

Product Code	CSB-PA015459LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43776
Immunogen	Recombinant Human Asparagine--tRNA ligase, cytoplasmic protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Asparagine--tRNA ligase, cytoplasmic (EC 6.1.1.22) (AsparaginyI-tRNA synthetase) (AsnRS), NARS
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NARS