



NARS Antibody, Biotin conjugated

Product Code	CSB-PA015459LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43776
Immunogen	Recombinant Human Asparagine--tRNA ligase, cytoplasmic protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Asparagine--tRNA ligase, cytoplasmic (EC 6.1.1.22) (AsparaginyI-tRNA synthetase) (AsnRS), NARS
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NARS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.