



NDUFV3 Antibody, HRP conjugated

Product Code	CSB-PA015671LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56181
Immunogen	Recombinant Human NADH dehydrogenase [ubiquinone] flavoprotein 3, mitochondrial protein (35-108AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NADH dehydrogenase [ubiquinone] flavoprotein 3, mitochondrial (Complex I-9kD) (CI-9kD) (NADH-ubiquinone oxidoreductase 9 kDa subunit) (Renal carcinoma antigen NY-REN-4), NDUFV3
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	NDUFV3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.