




NME5 Antibody

Product Code	CSB-PA015890LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56597
Immunogen	Recombinant Human Nucleoside diphosphate kinase homolog 5 protein (1-212AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoside diphosphate kinase homolog 5 (NDK-H 5) (NDP kinase homolog 5) (Inhibitor of p53-induced apoptosis-beta) (IPIA-beta) (Testis-specific nm23 homolog) (nm23-H5), NME5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NME5

Image



Immunofluorescent analysis of HeLa cells using CSB-PA015890LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.