



NME5 Antibody, HRP conjugated

Product Code	CSB-PA015890LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56597
Immunogen	Recombinant Human Nucleoside diphosphate kinase homolog 5 protein (1-212AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoside diphosphate kinase homolog 5 (NDK-H 5) (NDP kinase homolog 5) (Inhibitor of p53-induced apoptosis-beta) (IPIA-beta) (Testis-specific nm23 homolog) (nm23-H5), NME5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NME5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.