



NOX4 Antibody, HRP conjugated

Product Code	CSB-PA015961LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPH5
Immunogen	Recombinant Human NADPH oxidase 4 protein (210-424AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NADPH oxidase 4 (EC 1.6.3.-) (Kidney oxidase-1) (KOX-1) (Kidney superoxide-producing NADPH oxidase) (Renal NAD(P)H-oxidase), NOX4, RENOX
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	NOX4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.