



NPBWR2 Antibody, HRP conjugated

Product Code	CSB-PA015973LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48146
Immunogen	Recombinant Human Neuropeptides B/W receptor type 2 protein (1-45AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuropeptides B/W receptor type 2 (G-protein coupled receptor 8), NPBWR2, GPR8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NPBWR2