



Npr3 Antibody, Biotin conjugated

Product Code	CSB-PA016025LD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41740
Immunogen	Recombinant Rat Atrial natriuretic peptide receptor 3 protein (41-477AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Atrial natriuretic peptide receptor 3 (Atrial natriuretic peptide clearance receptor) (Atrial natriuretic peptide receptor type C) (ANP-C) (ANPR-C) (NPR-C), Npr3, Npr-c
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Npr3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.