



NR2C2AP Antibody, Biotin conjugated

Product Code	CSB-PA016053LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86WQ0
Immunogen	Recombinant Human Nuclear receptor 2C2-associated protein (7-139AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear receptor 2C2-associated protein (TR4 orphan receptor-associated 16 kDa protein), NR2C2AP, TRA16
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NR2C2AP