



NUDT14 Antibody, HRP conjugated

Product Code	CSB-PA016159LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95848
Immunogen	Recombinant Human Uridine diphosphate glucose pyrophosphatase protein (1-222AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uridine diphosphate glucose pyrophosphatase (UDPG pyrophosphatase) (UGPPase) (EC 3.6.1.45) (Nucleoside diphosphate-linked moiety X motif 14) (Nudix motif 14), NUDT14, UGPP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NUDT14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.