



NUDT21 Antibody, Biotin conjugated

Product Code	CSB-PA016170LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43809
Immunogen	Recombinant Human Cleavage and polyadenylation specificity factor subunit 5 protein (2-227AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cleavage and polyadenylation specificity factor subunit 5 (Cleavage and polyadenylation specificity factor 25 kDa subunit) (CFIm25) (CPSF 25 kDa subunit) (Nucleoside diphosphate-linked moiety X motif 21) (Nudix motif 21) (Pre-mRNA cleavage factor Im 25 kDa subunit), NUDT21, CFIM25 CPSF25 CPSF5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NUDT21
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.