



OIP5 Antibody, HRP conjugated

Product Code	CSB-PA016316LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43482
Immunogen	Recombinant Human Protein Mis18-beta protein (1-229AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein Mis18-beta (Cancer/testis antigen 86) (CT86) (Opa-interacting protein 5) (OIP-5), OIP5, MIS18B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	OIP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.