



OR1D2 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA016490LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P34982 |
| Immunogen | Recombinant Human Olfactory receptor 1D2 protein (159-196AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Olfactory receptor 1D2 (Olfactory receptor 17-4) (OR17-4) (Olfactory receptor OR17-6) (Olfactory receptor-like protein HGMP07E), OR1D2, OLF1R1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | OR1D2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |