



OR1D2 Antibody, Biotin conjugated

Product Code	CSB-PA016490LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P34982
Immunogen	Recombinant Human Olfactory receptor 1D2 protein (159-196AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactory receptor 1D2 (Olfactory receptor 17-4) (OR17-4) (Olfactory receptor OR17-6) (Olfactory receptor-like protein HGMP07E), OR1D2, OLFR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	OR1D2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.