



OR3A1 Antibody, FITC conjugated

Product Code	CSB-PA016636LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47881
Immunogen	Recombinant Human Olfactory receptor 3A1 protein (162-199AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactory receptor 3A1 (Olfactory receptor 17-40) (OR17-40) (Olfactory receptor OR17-15), OR3A1, OLFRA03
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	OR3A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.