



PAH Antibody, Biotin conjugated

Product Code	CSB-PA017396LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00439
Immunogen	Recombinant Human Phenylalanine-4-hydroxylase protein (2-452AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phenylalanine-4-hydroxylase (PAH) (EC 1.14.16.1) (Phe-4-monooxygenase), PAH
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PAH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.