



PCBD1 Antibody, HRP conjugated

Product Code	CSB-PA017514LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61457
Immunogen	Recombinant Human Pterin-4-alpha-carbinolamine dehydratase protein (2-104AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pterin-4-alpha-carbinolamine dehydratase (PHS) (EC 4.2.1.96) (4-alpha-hydroxy-tetrahydropterin dehydratase) (Dimerization cofactor of hepatocyte nuclear factor 1-alpha) (DCoH) (Dimerization cofactor of HNF1) (Phenylalanine hydroxylase-stimulating protein) (Pterin carbinolamine dehydratase) (PCD), PCBD1, DCOH PCBD
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PCBD1