



PCBD1 Antibody, Biotin conjugated

Product Code	CSB-PA017514LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61457
Immunogen	Recombinant Human Pterin-4-alpha-carbinolamine dehydratase protein (2-104AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pterin-4-alpha-carbinolamine dehydratase (PHS) (EC 4.2.1.96) (4-alpha-hydroxy-tetrahydropterin dehydratase) (Dimerization cofactor of hepatocyte nuclear factor 1-alpha) (DCoH) (Dimerization cofactor of HNF1) (Phenylalanine hydroxylase-stimulating protein) (Pterin carbinolamine dehydratase) (PCD), PCBD1, DCOH PCBD
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PCBD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.