



# PCBD1 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA017514LD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P61457
<b>Immunogen</b>	Recombinant Human Pterin-4-alpha-carbinolamine dehydratase protein (2-104AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Pterin-4-alpha-carbinolamine dehydratase (PHS) (EC 4.2.1.96) (4-alpha-hydroxy-tetrahydropterin dehydratase) (Dimerization cofactor of hepatocyte nuclear factor 1-alpha) (DCoH) (Dimerization cofactor of HNF1) (Phenylalanine hydroxylase-stimulating protein) (Pterin carbinolamine dehydratase) (PCD), PCBD1, DCOH PCBD
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	PCBD1