



PCP4 Antibody

Product Code	CSB-PA017636HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48539
Immunogen	Recombinant Human Calmodulin regulator protein PCP4 protein (1-62AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calmodulin regulator protein PCP4 (Brain-specific polypeptide PEP-19) (Purkinje cell protein 4), PCP4, PEP19
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PCP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.