



Pcsk5 Antibody, FITC conjugated

Product Code	CSB-PA017644LC01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q04592
Immunogen	Recombinant Mouse Proprotein convertase subtilisin/kexin type 5 protein (117-452AA)
Raised In	Rabbit
Species Reactivity	Mouse
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proprotein convertase subtilisin/kexin type 5 (EC 3.4.21.-) (Proprotein convertase 5) (PC5) (Proprotein convertase 6) (PC6) (Subtilisin-like proprotein convertase 6) (SPC6) (Subtilisin/kexin-like protease PC5), Pcsk5
Immunogen Species	Mus musculus (Mouse)
Research Area	Epigenetics and Nuclear Signaling
Target Names	Pcsk5