



Pcsk5 Antibody, Biotin conjugated

Product Code	CSB-PA017644LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q04592
Immunogen	Recombinant Mouse Proprotein convertase subtilisin/kexin type 5 protein (117-452AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proprotein convertase subtilisin/kexin type 5 (EC 3.4.21.-) (Proprotein convertase 5) (PC5) (Proprotein convertase 6) (PC6) (Subtilisin-like proprotein convertase 6) (SPC6) (Subtilisin/kexin-like protease PC5), Pcsk5
Immunogen Species	Mus musculus (Mouse)
Research Area	Epigenetics and Nuclear Signaling
Target Names	Pcsk5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.