



PCTP Antibody, Biotin conjugated

Product Code	CSB-PA017651LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKL6
Immunogen	Recombinant Human Phosphatidylcholine transfer protein (1-214AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylcholine transfer protein (PC-TP) (START domain-containing protein 2) (StARD2) (StAR-related lipid transfer protein 2), PCTP, STARD2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PCTP