



PDCD1LG2 Antibody, Biotin conjugated

Product Code	CSB-PA017667ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BQ51
Immunogen	Recombinant Human Programmed cell death 1 ligand 2 protein (20-220AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Programmed cell death 1 ligand 2 (PD-1 ligand 2) (PD-L2) (PDCD1 ligand 2) (Programmed death ligand 2) (Butyrophilin B7-DC) (B7-DC) (CD antigen CD273), PDCD1LG2, B7DC CD273 PDCD1L2 PDL2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PDCD1LG2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.