



PDIA4 Antibody, Biotin conjugated

Product Code	CSB-PA017722LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13667
Immunogen	Recombinant Human Protein disulfide-isomerase A4 protein (101-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein disulfide-isomerase A4 (EC 5.3.4.1) (Endoplasmic reticulum resident protein 70) (ER protein 70) (ERp70) (Endoplasmic reticulum resident protein 72) (ER protein 72) (ERp-72) (ERp72), PDIA4, ERP70 ERP72
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDIA4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.