



PGLYRP4 Antibody, FITC conjugated

Product Code	CSB-PA017865LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96LB8
Immunogen	Recombinant Human Peptidoglycan recognition protein 4 protein (18-317AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidoglycan recognition protein 4 (Peptidoglycan recognition protein I-beta) (PGLYRP4) (PGRP-I-beta) (Peptidoglycan recognition protein intermediate beta), PGLYRP4, PGRPIB
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PGLYRP4