



PGLYRP4 Antibody, FITC conjugated

Product Code	CSB-PA017865LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96LB8
Immunogen	Recombinant Human Peptidoglycan recognition protein 4 protein (18-317AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidoglycan recognition protein 4 (Peptidoglycan recognition protein I-beta) (PGLYRPIbeta) (PGRP-I-beta) (Peptidoglycan recognition protein intermediate beta), PGLYRP4, PGRPIB
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PGLYRP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.