



PGLYRP4 Antibody, Biotin conjugated

| | |
|---------------------|---|
| Product Code | CSB-PA017865LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96LB8 |
| Immunogen | Recombinant Human Peptidoglycan recognition protein 4 protein (18-317AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Peptidoglycan recognition protein 4 (Peptidoglycan recognition protein I-beta) (PGLYRP4) (PGRP-I-beta) (Peptidoglycan recognition protein intermediate beta), PGLYRP4, PGRPIB |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | PGLYRP4 |