



PKLR Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA018071LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P30613 |
| Immunogen | Recombinant Human Pyruvate kinase PKLR protein (1-574AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Pyruvate kinase PKLR (EC 2.7.1.40) (Pyruvate kinase 1) (Pyruvate kinase isozymes L/R) (R-type/L-type pyruvate kinase) (Red cell/liver pyruvate kinase), PKLR, PK1 PKL |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | PKLR |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |