



PKLR Antibody, Biotin conjugated

Product Code	CSB-PA018071LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30613
Immunogen	Recombinant Human Pyruvate kinase PKLR protein (1-574AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pyruvate kinase PKLR (EC 2.7.1.40) (Pyruvate kinase 1) (Pyruvate kinase isozymes L/R) (R-type/L-type pyruvate kinase) (Red cell/liver pyruvate kinase), PKLR, PK1 PKL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PKLR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.