



PLK1 Antibody, HRP conjugated

Product Code	CSB-PA018193LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53350
Immunogen	Recombinant Human Serine/threonine-protein kinase PLK1 protein (1-603AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase PLK1 (EC 2.7.11.21) (Polo-like kinase 1) (PLK-1) (Serine/threonine-protein kinase 13) (STPK13), PLK1, PLK
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PLK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.