



PLK1 Antibody, Biotin conjugated

Product Code	CSB-PA018193LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53350
Immunogen	Recombinant Human Serine/threonine-protein kinase PLK1 protein (1-603AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase PLK1 (EC 2.7.11.21) (Polo-like kinase 1) (PLK-1) (Serine/threonine-protein kinase 13) (STPK13), PLK1, PLK
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PLK1