




PLOD2 Antibody

Product Code	CSB-PA018200LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00469
Immunogen	Recombinant Human Procollagen-lysine,2-oxoglutarate 5-dioxygenase 2 protein (301-600AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Procollagen-lysine,2-oxoglutarate 5-dioxygenase 2 (EC 1.14.11.4) (Lysyl hydroxylase 2) (LH2), PLOD2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PLOD2

Image



Western Blot

All lanes: PLOD2 antibody at 12µg/ml

Lane 1: U87 whole cell lysate

Lane 2: U251 whole cell lysate

Lane 3: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 85, 88, 50 kDa

Observed band size: 85 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.