



PNLIPRP3 Antibody, Biotin conjugated

Product Code	CSB-PA018265LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q17RR3
Immunogen	Recombinant Human Pancreatic lipase-related protein 3 protein (201-467AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pancreatic lipase-related protein 3 (PL-RP3) (EC 3.1.1.3), PNLIPRP3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PNLIPRP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.