



POGK Antibody, FITC conjugated

Product Code	CSB-PA018298LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P215
Immunogen	Recombinant Human Pogo transposable element with KRAB domain protein (101-400AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pogo transposable element with KRAB domain, POGK, KIAA1513
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	POGK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.