



POLB Antibody, HRP conjugated

Product Code	CSB-PA018302LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06746
Immunogen	Recombinant Human DNA polymerase beta protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA polymerase beta (EC 2.7.7.7) (EC 4.2.99.-), POLB
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	POLB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.