



# POLR2K Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA018340LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P53803
<b>Immunogen</b>	Recombinant Human DNA-directed RNA polymerases I, II, and III subunit RPABC4 protein (1-58AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	DNA-directed RNA polymerases I, II, and III subunit RPABC4 (RNA polymerases I, II, and III subunit ABC4) (ABC10-alpha) (DNA-directed RNA polymerase II subunit K) (RNA polymerase II 7.0 kDa subunit) (RPB7.0) (RPB10alpha), POLR2K
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	POLR2K
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.