



POLR2K Antibody, Biotin conjugated

Product Code	CSB-PA018340LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53803
Immunogen	Recombinant Human DNA-directed RNA polymerases I, II, and III subunit RPABC4 protein (1-58AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-directed RNA polymerases I, II, and III subunit RPABC4 (RNA polymerases I, II, and III subunit ABC4) (ABC10-alpha) (DNA-directed RNA polymerase II subunit K) (RNA polymerase II 7.0 kDa subunit) (RPB7.0) (RPB10alpha), POLR2K
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	POLR2K
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.