



POMP Antibody, Biotin conjugated

Product Code	CSB-PA018365LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y244
Immunogen	Recombinant Human Proteasome maturation protein (1-141AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteasome maturation protein (Proteasemblin) (Protein UMP1 homolog) (hUMP1) (Voltage-gated K channel beta subunit 4.1), POMP, C13orf12 UMP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	POMP